

## **KRAS Somatic Mutation Analysis, Peritoneal Fluid**

**Test ID:** KRASW

**Explanation:** On the effective date, the specimen stability will be updated.

<b>Current Specimen Stability</b>			
<b>Specimen</b>	<b>Temperature</b>	<b>Time</b>	<b>Special Container</b>
Varies	Refrigerated		

<b>New Specimen Stability</b>			
<b>Specimen</b>	<b>Temperature</b>	<b>Time</b>	<b>Special Container</b>
Peritoneal	Refrigerated (preferred)	10 days	
	Ambient	5 days	

### **Questions**

Contact Michelle Rath, Laboratory Resource Coordinator at 800-533-1710.