
MVista Blastomyces Quantitative Antigen, Fluid

Test ID: FBMO

Explanation: Effective February 10, 2025, test ID FBMO, referred to Miravista Diagnostics, will reflect the following changes:

Current Cautions
Cross-reactions are seen with histoplasmosis, paracoccidioidomycosis, penicilliosis, less frequently in coccidioidomycosis, rarely in aspergillosis and possibly sporotrichosis. Sputolysin, sodium hydroxide, and potassium hydroxide treatment degrade the analyte detected in the assay

New Cautions
When tested in cultures of 10^5 (1,000,000)– 10^6 (10,000,000) organisms/mL, cross-reactions occurred with <i>Histoplasma</i> spp., <i>Coccidioides</i> spp., <i>Paracoccidioides brasiliensis</i> , <i>Talaromyces marneffe</i> , <i>Aspergillus nidulans</i> , and <i>Candida tropicalis</i> .

Current NY State Information
No

New NY State Information
Yes

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.