
Risperidone (Risperdal) and 9-Hydroxyrisperidone

Test ID: FRISP

Explanation:

Effective immediately, testing is delayed for test ID: FRISP, referred to Medtox Laboratories, due to QC issues. There is no ETA currently. Medtox will stabilize samples and test on a first in, first out basis when able.

Specimen Stability Information:

Specimen Type	Temperature	Time
Varies	Refrigerated (preferred)	7 days
	Ambient	72 hours
	Frozen	180 days

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.