

Hemoglobin Variant, A2 and F Quantitation, Blood

Test ID: HGBCE

Explanation: On the effective date, sodium heparin (green top) tubes will no longer be acceptable.

Current Specimen Required
<p>Preferred: Lavender top (EDTA) Acceptable: Yellow top (ACD), green top (heparin) Specimen Volume: 4 mL Collection Instructions: 1. Submit fresh specimen. 2. Send whole blood specimen in original tube. Do not aliquot.</p>

New Specimen Required
<p>Preferred: Lavender top (EDTA) Acceptable: Yellow top (ACD) Specimen Volume: 4 mL Collection Instructions: 1. Submit fresh specimen. 2. Send whole blood specimen in original tube. Do not aliquot.</p>

Questions

Contact Melissa Lonzo, Laboratory Resource Coordinator at 800-533-1710.