

Influenza Virus Type A and Type B, and Respiratory Syncytial Virus (RSV), Molecular Detection, PCR, Varies

Test ID: FLUMS; performed at Mayo Clinic Laboratories Florida.

Explanation:

In an effort to realign current respiratory virus test offerings, this test will be made obsolete on the effective date.

Recommended Alternative Test:

Severe Acute Respiratory Syndrome Coronavirus 2 (SARS CoV-2), Influenza A and B, and Respiratory Syncytial Virus, Molecular Detection, Bronchoalveolar Lavage

Test ID: SFRB; performed at Mayo Clinic Laboratories Florida.

Methodology:

Real-Time Reverse Transcription Polymerase Chain Reaction (RT-PCR)

Reference Values:

Undetected

Specimen Requirements:

Source:	Bronchoalveolar lavage (BAL)
Container/Tube:	Sterile container
Specimen Volume:	1 mL
Additional Information:	Do not aliquot into viral transport media.

Specimen Stability Information:

Specimen Type	Temperature	Time
Bronchoalveolar Lavage	Refrigerated (preferred)	72 hours
	Ambient	24 hours

CPT Code:

0241U

Day(s) Performed: Monday through Sunday **Report Available:** Same day/1 to 3 days**Questions**

Contact Bonnie Meyers, Laboratory Resource Coordinator at 800-533-1710.