

Pyridostigmine, Serum/Plasma

Test ID: FPYD

Explanation: Effective May 19, 2025, test ID FPYD, referred to NMS Labs, will become obsolete.

Recommended Alternative Test:

Miscellaneous NMS Testing

Test ID: ZW86

Referral Lab Code: 4038SP

Methodology:

Gas Chromatography (GC)

Reference Values:

None Detected

Specimen Requirements:

Submit only 1 of the following specimens:

Specimen Type: Serum

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 5 mL

Collection Instructions:

1. Draw blood in a plain, red-top tube(s). Serum gel tube is not acceptable.
2. Centrifuge and freeze immediately. Send 5 mL of serum frozen in a plastic, preservative-free vial.
3. Label specimen appropriately (serum).

Specimen Type: Plasma

Collection Container/Tube:

Preferred: Lavender-top EDTA)

Acceptable: Green-top (Heparin)

Submission Container/Tube: Plastic vial

Specimen Volume: 5 mL

Collection Instructions:

1. Draw blood in a lavender-top (EDTA) or green-top (Heparin) tube(s). Plasma gel tube is not acceptable.
2. Centrifuge and freeze immediately. Send 5 mL of plasma frozen in a preservative-free plastic vial.
3. Label specimen appropriately (plasma).

Minimum volume: 2.2 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Varies	Frozen	21 Days

CPT Code:

80375

Fee:

\$404.00

Day(s) Setup: Monday-Sunday

Analytic Time: 9-11 days

Questions

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.