



Leukotriene E4, 24 Hour, Urine

Test ID: TLTE4

Explanation:

Updating specimen stability and adding a frozen creatinine code to match the assay's preferred temperature.

Current Profile Information	
Test ID	Reporting Name
LTE4T	Leukotriene E4, 24 Hr, U
CRT24	Creatinine, 24 HR, U

New Profile Information	
Test ID	Reporting Name
LTE4T	Leukotriene E4, 24 Hr, U
CRT2F	Creatinine, 24 HR, U

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Frozen (Preferred)	30 days
	Ambient	24 hours
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Frozen (Preferred)	28 days
	Ambient	24 hours
	Refrigerated	7 days

Questions

Contact Joshua Yang, Laboratory Resource Coordinator at 800-533-1710.