

Immunofixation, CSF

Test ID: FIXCF

Explanation: Per notification received from Quest Diagnostics Chantilly, Test ID FIXCF, will experience the following changes, effective immediately.

Current Specimen Stability		
Specimen	Temperature	Time
CSF	Refrigerated (preferred)	7 days
	Ambient	48 hours
	Frozen	180 days

New Specimen Stability		
Specimen	Temperature	Time
CSF	Frozen (preferred)	30 days
	Ambient	6 days
	Refrigerated	10 days

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.