

Alzheimer Disease Evaluation, Spinal Fluid

Test ID: ADEVL

Explanation: Specimen requirement and reference changes as reflected below.

Current Reference Value	
Beta-amyloid (1-42) (Abeta42):	>1026 pg/mL
Total-Tau:	< or = 238 pg/mL
Phosphorylated-Tau 181:	< or = 21.7 pg/mL
p-Tau/Abeta42:	< or = 0.023

New Reference Value	
Beta-amyloid (1-42) (Abeta42):	>834 pg/mL
Total-Tau:	< or = 238 pg/mL
Phosphorylated-Tau 181:	< or = 21.6 pg/mL
p-Tau/Abeta42:	< or = 0.028

Current Specimen Required	
Container/Tube:	
Preferred:	CSF AD Biomarker Tube (2.5 mL)
Acceptable:	Sarstedt 72.703.600 (1.5 mL) or Sarstedt 72.694.600 (2 mL)

New Specimen Required	
Container/Tube:	
Preferred:	Sarstedt CSF False Bottom Tube 63.614.625 (2.5 mL)
Acceptable:	Sarstedt 72.703.600 (1.5 mL) or Sarstedt 72.694.600 (2 mL)

Current Requested Volume
1.5 to 2 mL

New Requested Volume
1.5 to 2.5 mL

Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.