



***Entamoeba histolytica* Antibody, Serum**

Test ID: EHOLG

Explanation:

Specimen stability for *Entamoeba histolytica* Antibody, Serum will be updated as below.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	30 days
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	30 days

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.