



Supersaturation Profile, Self-Collect, 24 Hour, Urine

Test ID: SSH24

Explanation: On the effective date, ambient and frozen samples will be accepted within the stability listed below, and the Days and Times Performed will be updated.

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated	14 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated	14 days
	Frozen	14 days
	Ambient	3 days

Current Days Performed
Monday through Friday

New Days Performed
Monday through Sunday

Current Times Performed
8 a.m. – 4 p.m.

New Times Performed
Monday through Friday; 7 a.m.-4 p.m., Saturday and Sunday 8 a.m.- 3:30 p.m

Questions

Contact Nancy Benson, Laboratory Resource Coordinator at 800-533-1710.